

mbo_utilities

Cheat Sheet - Python API, CLI & GUI

Python API: Core I/O

`imread(path)`
Lazy-load any supported format

`imwrite(arr, path, ext)`
Stream-write to disk

`get_metadata(path)`
Extract file metadata dict

`get_voxel_size(path)`
Get physical dimensions (μm)

`get_files(path, **kw)`
Discover files with filtering

Examples

```
# Load data
from mbo_utilities import imread
arr = imread('/path/to/data.tiff')
arr = imread('/path/to/raw/') # ScanImage
```

```
# Write data
from mbo_utilities import imwrite
imwrite(arr, 'out.zarr', ext='.zarr')
imwrite(arr, 'out/', planes=[0,1,2])
```

```
# Get metadata
from mbo_utilities import get_metadata
meta = get_metadata('/path/to/data')
print(meta['nframes'], meta['shape'])
```

Utilities & Visualization

`save_mp4(fname, images, **kw)`
Export video from 3D array

`save_png(fname, data)`
Save image via matplotlib

`norm_minmax(images)`
Normalize to 0-1 range

`smooth_data(data, window)`
Temporal smoothing

`subsample_array(arr, factor)`
Downsample array

`files_to_dask(files)`
Build Dask array from files

`expand_paths(paths)`
Expand wildcards/lists

Examples

```
# Video export
from mbo_utilities import save_mp4
save_mp4('movie.mp4', arr[:500],
        framerate=30, temporal_avg=5)
```

```
# Dask arrays
from mbo_utilities import files_to_dask
darr = files_to_dask(tiff_files,
                    chunk_t=250)
```

CLI Commands

<code>mbo view [PATH]</code>	Launch GUI viewer	<code>--roi 0,1 --widget --metadata</code>
<code>mbo convert IN OUT</code>	Convert formats	<code>-e .zarr -p 0,1,2 --fix-phase --register-z</code>
<code>mbo info PATH</code>	Show file metadata	<code>--metadata</code>
<code>mbo scanphase [PATH]</code>	Analyze scan phase	<code>-o output/ --format png --show</code>
<code>mbo formats</code>	List supported formats	
<code>mbo --download-notebook</code>	Get user guide	<code>[PATH]</code>

Quick Examples

```
mbo /path/to/data.tiff          # View TIFF in GUI
mbo convert raw/ out.zarr -e .zarr # Convert to Zarr
mbo convert data.tiff out/ --fix-phase # Fix bidirectional scan
mbo view data/ --roi 0,1         # View specific ROIs
```

Supported Formats

Input

<code>.tif, .tiff</code>	BigTIFF, OME-TIFF, ScanImage
<code>.zarr</code>	Zarr v3, OME-NGFF
<code>.h5, .hdf5</code>	HDF5 datasets
<code>.bin</code>	Suite2p binary + ops.npy
<code>.npy</code>	NumPy arrays
<code>In-memory</code>	NumPy/Dask arrays

Output

<code>.tiff</code>	BigTIFF (streaming write)
<code>.zarr</code>	Zarr v3 with OME metadata
<code>.h5</code>	HDF5 with chunking
<code>.bin</code>	Suite2p binary format
<code>.npy</code>	NumPy array
<code>.mp4</code>	Video export

Lazy Array Types (returned by imread)

- `MboRawArray` - Raw ScanImage multi-ROI data with metadata
- `TiffArray` - Memory-mapped TIFF access
- `ZarrArray` - Chunked cloud-ready arrays
- `Suite2pArray` - Suite2p binary with ops integration

GUI Features

Preview & Visualization

- Image Viewer - FastPlotLib 2D/3D rendering with WGPU
- Frame Navigation - Time slider with playback controls
- Z-Plane Slider - Navigate through imaging planes
- Window Functions - mean, max, std, mean-subtracted
- Scan-Phase Correction - Fix bidirectional artifacts
- Contrast Controls - V-Min/V-Max adjustment
- Summary Stats - Per-plane mean, std, SNR tables

Processing & Export

- Spatial Crop - Select ROI region for processing
- Suite2p Pipeline - Integrated registration & cell detection
- Registration Settings - Rigid/non-rigid, 1P mode options
- Save As Dialog - Export to .tiff/.zarr/.h5/.bin
- Multi-ROI Support - Process ROIs separately or combined
- Suite3D Registration - Axial z-plane alignment
- Phase Correction - Automatic bidirectional scan fix

ROI Diagnostics (Suite2p Results)

- dF/F Traces - Adjustable baseline (median/percentile)
- Quality Metrics - SNR, skewness, activity histograms
- Filter Sliders - Interactive threshold adjustment
- Auto-save - Syncs iscell.npy to disk on change
- File Watching - Detects external modifications
- Suite2p Sync - Bi-directional with Suite2p GUI
- ROI Statistics - Detailed per-ROI information

Launching the GUI

From command line

```
mbo view /path/to/data  
mbo /path/to/data.tiff
```

From Python

```
from mbo_utilities import run_gui  
run_gui('/path/to/data')  
run_gui() # opens file dialog
```